

## Automated Flow Chemisorption Analyzer request form for External User

**Instrument Name:** Chemisorption analyser (Anton Paar) Model: - ChemBET Pulsar

---

### **Instructions:**

1. **Sample Form:** Powder (~100 mg).
  2. **Degassing Temperature Range:** From Room Temperature to 300 °C. The analysis cannot be carried out if degassing and analysis conditions are not specified.
  3. TGA data is required
  4. **Prohibited Samples:** Explosive, poisonous samples, or those releasing toxic/obnoxious gases/fumes upon heating cannot be analysed.
  5. The user must ensure that samples will not melt or change during degassing.
  6. **Repeat Analysis:** A new form must be filled out for each repeat analysis, and the sample will be analysed as per the queue.
- 

### **Personal Information**

1. **Type of Industry:** Academic / Industrial / National Lab
  2. **Name of the User:** \_\_\_\_\_
  3. **User Email ID:** \_\_\_\_\_
  4. **Registration No:** \_\_\_\_\_
  5. **Mobile Number:** \_\_\_\_\_
  6. **Department / Institution name:** \_\_\_\_\_
  7. **GST Number:** \_\_\_\_\_
  8. **Number of Samples (Only 1 at a time)** \_\_\_\_\_
  9. **Sample Code:** \_\_\_\_\_
  10. **Degassing Temperature.** \_\_\_\_\_
  11. **Degassing duration (in hours):** \_\_\_\_\_
  12. **Type of Analysis:** TPR/TPD/TPO
  13. **Additional Information (If any):**
- 

**Note: Final data will be shared in ASCII and PDF format to your email.**