INDIAN INSTITUTE OF TECHNOLOGY, BOMBAY

Sophisticated Analytical Instrument Facility (SAIF)

Environmental Scanning Electron Microscope (ESEM)

Analysis Request Form

Applicant Details	Applicant Details		
User belongs to:			
User name:			
Institute/University/Organisation:			
Email ID: Mobile No.:			
Name of Guide/PI:			
Guide/PI Email ID: Guide/PI Mobile No.:			
Address of Institute/Organization:			
Sample information:			
Number of samples			
Sample code			
Sample type	Biological / Composite Material / Thin Film / Crystalline Solid / Metal/Polymer/Ceramic/Composite/ Other (Please specify):* *biological samples should be submitted after primary fixation		
Detailed description of the sample (Refer Annexure I before filling)			
Sample form	Powder/Pellet/Suspension/On glass substrate/Film/Others		
If sample is Pellet, thin film (Refer Annexure II before filling)	Mention sample dimensions		
If sample is powder, Dispersion required (Refer Annexure II before filling)	Yes/No If Yes, Medium for dispersion		
Nature of the sample	Non-Conducting / Conducting / Semi-Conductor		
Sample is	Magnetic /Non Magnetic		
Type of analysis (Kindly tick):			
SEM Image only	SEM Image + EDS	Only EDS	

SEM Analysis requirement:

SEM Image	Secondary Electron Image / Back Scattered Electron Image
Sample to be mounted	Planar/Cross Section/ Powder directly loaded on sample holder/ Drop cast for liquid or suspension samples
Expected Morphology	Brief description of shape
Expected Particle Size	
If EDS is required	Point EDS /area EDS/EDS mapping/ line scan List of expected elements (for EDS):

Material safety data:

If you are submitting more than one sample which are different in nature/composition, submit separate MSDS

Sample Properties	Carcinogenic (level) Toxic Radioactive Corrosive Explosive Flammable Other (specify):
Moisture	Present Absent NA
Volatile organic compound	Present Absent NA
Stability of sample	Stable under RTP Hygroscopic Sublimes Reactive in: Air Light Heat Vacuum Moisture May decompose when exposed to accelerated electron beam
Mention the storage and handling conditions if anything specific	
Whether incompatible with any material-	Yes No Specify the materials):
Health hazards	Yes No (irritant to skin/irritant to eyes/harmful to skin/toxic if inhaled/toxic if ingested)
First aid measures	Eye/Skin/Inhalation/ Ingestion/Others (specify):
Disposal Method of sample	
Please fill appropriate numbers in the NFPA diamond: (*reference image attached below)	
Additional information if any	

*Along with this form MSDS should be submitted if available.

Note: All Samples will be discarded after 15 days of analysis. If you wish to collect the samples then you are required to make arrangement for the same. SAIF office will not dispatch the same to users under any circumstances

Declaration

I confirm that the samples submitted for analysis are for research purpose only and the above furnished details are correct and true to the best of my knowledge. I understand that I will be held responsible for any damages arising from incorrect information provided by me against material safety data.

I agree to acknowledge DST and SAIF/CRNTS, IIT Bombay for providing (Instrument name) analytical facility for my research work, in my publications. I also agree to send the publication reference (Journal name/volume number/names of the authors/date of issue of the publication etc) to office.saif@iitb.ac.in

I declare that the "Content of this report is meant for our information only and we will not use the content of this report for advertisement, evidence, litigation or quote as certificate to third party" I accept that all the issued reports/results (Soft/hard) will not carry any Signature or Seal and Stamp of SAIF/CRNTS IIT Bombay.

Signature of the User

Signature of the In Charge/HOD/PI with College/P.I. Guide seal / stamp

Date: Place:

* Reference image for filling NFPA diamond:



Annexure I, II is for your reference (kindly do not print)

Annexure I

For filling detailed description of the sample: kindly refer to the below sub categories and examples. If your sample details do not match with the below list, please give the correct sample type and sample description.

Sample type	Description
Biological	Cells (Give type), Tooth, Gels, Scaffolds, Bone, Biofilm, Tissue,
	Leaf/plant extracts, Insect/Insect parts, Lipids/Liposome's,
	Proteins, Blood cells, Bacteria (Give type), Sludge, Fibrin gel,
Polymer	Resin, Alginate, Polystyrene, Polypropylene, PDMS, PVC,
	Polymeric microspheres, Fibers, Thermoplastic polyurethanes,
	Polymeric scaffold
Metal	Alloy, Chips, Micro tools, Fractography
Geological	Soil, Flyash, Sand, Activated carbon, Brick, Cement
Nanomaterials	CNT, Nano particles (Give type), Ferrite
Thin film	Specify the material:, Substrate: Glass/ Copper/Conducting
	material/Silicon wafer
Ceramic or Composite material	Detailed description of the sample/composite material

EXAMPLES:

Sample type	Description
Biological	Shrimp waste extract
Biological and Nanomaterials	Au/Ag nano particles prepared from plant extract
Biological	Cells, blood cells or animal cells or E.coli/Staph. aureus
Polymer	polystyrene nanoparticles
Thinfilm	Material: ZnO/TiO2/CZTS substrate: glass
Nanomaterials	Gold nano particles or CNT or
Composite material	CNT in polymer, CNT/ carbon and graphene

^{**}For Powder sample: Mentioning sample type and sample description is mandatory.

Annexure II

Sample preparation instructions for SEM:

- **Medium for Dispersion** Ethanol / Methanol / Water /Iso-propyl alcohol. Any other medium should be provided by the user. Dispersion will be done by ultrasonication.
- Kindly mark the edge of the sample to be observed for Cross section.
- Base of the sample should be flat for mounting on sample holder.
- Biological samples will be accepted only after primarily fixation with suitable fixative.
- Sample preparation if any should be done at user end (cutting the sample for CS, freeze fracturing, sample fixation for biological samples, staining of samples, oven drying should be done by the user)
- For any further query, kindly contact on Email: esem@iitb.ac.in, Contact: 022-2159-6862

IMPORTANT NOTE:

- Potentially hazardous/toxic/radioactive samples may not be accepted for analysis.
- Requisition letter addressing Head, SAIF/CRNTS along with completed form should be send by post or submitted in person to **SAIF/CRNTS office, IIT Bombay, Powai, Mumbai-400076.**
- Your appointment will be as per the queue, once we receive the Requisition letter, duly filled form and advance payment. Partially filled form will not be registered.
- We prefer that you/ your representative, who know/understands the sample/material and what is expected to be seen, will be present on the day of appointment. If the user is not present representative data will be taken for the samples.
- Attach reference images for the sample (if any) with the form.