

- **Appointment:** The users will be informed about their date and time-slot by email.
- **Sample Submission:** Samples are to be submitted at the time of registration or brought along on the date of your appointment for your sample analysis. Users can send the samples via courier also
- **Results:** After the slot is done the results will be sent by email for external users.

GENERAL INSTRUCTIONS TO EXTERNAL USERS

- ✚ The slot is of 2 hours and maximum 4 number of samples are allowed in a slot, the charges mentioned are for per slot.
- ✚ The users should visit <http://rnd.iitb.ac.in> website for the sample preparation instructions as well as for the contact details and location of the Facility
- ✚ The experimental data provided is only for research / development purposes.
These cannot be used as certificates in legal disputes.



INDIAN INSTITUTE OF TECHNOLOGY BOMBAY, POWAI, MUMBAI

IRCC Confocal Laser Scanning Microscopy Central Facility,
Department of Biosciences and Bioengineering
Email: lsmconfocal@iitb.ac.in Ph: 022-21596746

Registration Form for External Users

- 1) Name of the User:
2) Name of the Institute/Organization:
(Appropriate letter to be enclosed; please see instruction sheet)
3) GSTIN:
4) Email and Tel No:
5) Name of Research Supervisor:

Sample Information

- 6) Description of the sample to be imaged:
(E.g. fixed mammalian cells on a slide, live bacteria in a petridish, etc. Please give as many details as possible)
7) Number of samples:
8) Fluorophores present in the sample:
9) Type of measurement requested:
Multi color confocal Imaging with DIC Z-stack imaging Time series Multi-point Imaging
FRET FRAP FLIM FCS Tile Imaging
10) Special Requirement if any:

We agree to acknowledge "Confocal Laser Scanning Microscopy Central Facility of IRCC" of IIT BOMBAY in our Publications/Reports/Thesis in which the data is used with due feedback through email lsmconfocal@iitb.ac.in. Please send us a copy of the publication/thesis report.

Signature of the User:

Note: 1) Samples should be brought by the user at the allotted time and the user should be present while the samples are being examined.
2) Fixed samples should be sealed between glass side and cover-slip and for live samples it should be in 35 mm glass bottom petri-dishes.

For Official Use only

External User Registration No: Date:

Table with 2 columns and 2 rows: Date of experiment, Name of the operator, Results sent date, Signature of operator.