



IIT Bombay Nanofabrication Facility

Tool Name: Plasma Asher System (EUR-200)

Nano 2 Lab

Policy Document

Updated On: 28th October, 2024

General Policy:

- Use lint free cloth and IPA only for cleaning chamber.
- If you want to clean the instrument body, use an IPA soaked cloth.
- If any problem arises in the system while you are working:
- Please make a note of it in the logbook Call the system owner/email and inform him. Place a note near the system.
- Please do not try to rectify the problem, unless you are with an experienced AU or SO or you have been instructed by the SO to do so.

Training Policy:

- Minimum Two Eyes On session required before Hands On.
- Minimum Five Hands On session required before Authorization Test.
- Authorization Test will be conducted by System Owner or Assigned Authorized User.
- If User cleared the test would become an Authorized User.
- If fails then minimum Three Hands On session required before Re-Authorization Test.
- If a user's authorization gets revoked due to long inactivity then he has to give Authorization Test after completing three Hands On sessions with AUs

